

ADVANCED
CORE DECOMPRESSION SYSTEM

X-REAM[®]

Percutaneous
Expandable Reamer

PRO-DENSE[®]

Core Decompression
Procedure Kit

SURGICAL TECHNIQUE



 **WRIGHT**[™]
FOCUSED EXCELLENCE

Advanced Core Decompression System

Indications*

PRO-DENSE® resultant paste is intended for use as a bone graft substitute to be injected or digitally packed into open bone voids/gaps that are not intrinsic to the stability of bony structure of the skeletal system (i.e., the extremities and pelvis) to cure in situ. These open bone voids may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The paste provides a bone graft substitute that resorbs and is replaced with bone during the healing process.

The PRO-DENSE® paste cured in situ provides an open void/gap filler that can augment provisional hardware (e.g. KWires) to help support bone fragments during the surgical procedure. The cured paste acts only as a temporary support media and is not intended to provide structural support during the healing process.

PRO-DENSE® is provided sterile for single use only.

*US Indications

For the PRO-DENSE® Core Decompression Procedure Kit:

The PRO-DENSE® Core Decompression Procedure Kit, consisting of a bone void filler and manual surgical instruments, is intended to be used during core decompression procedures. The bone void filler component resorbs and is replaced with bone during the healing process. The bone void filler included in the PRO-DENSE® Core Decompression Procedure Kit is not intended to be used as a load-bearing device.

Contraindications

The PRO-DENSE® Bone Graft Substitute injectable paste is contraindicated where the device is intended as structural support in load-bearing bone and in articulating surfaces. Conditions representing relative contraindications include:

- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative bone disease
- Closed bone void/gap
- Pregnancy
- Uncooperative patients who will not or cannot follow postoperative instructions, including individuals who abuse drugs and/or alcohol
- Hypercalcemia
- Renal compromised patients
- Patients with a history of or active Pott's disease

X-REAM® Percutaneous Expandable Reamer

Reusable



PRO-DENSE® Core Decompression Procedure Kit



Advanced Core Decompression

Instrumentation and Grafting

The Advanced Core Decompression System includes the reusable X-REAM® Percutaneous Expandable Reamer that allows optimized debridement when used in conjunction with the PRO-DENSE® Core Decompression Procedure Kit (CDK) to prepare for a standard core decompression. The procedure kit, sold separately, includes single-use, disposable instruments that are designed to efficiently facilitate a standard core decompression, and PRO-DENSE® Injectable Graft for backfilling the surgically-created defect.

These instruments have been carefully selected and tested to simplify the technique for efficiency and consistency and to possibly provide a more cost effective outcome.¹

Note: *The PRO-DENSE® Core Decompression Procedure Kit is designed for single site usage.*

Minimally-Invasive Technique

When properly used, the expandable reamer tool allows optimal debridement of dead bone through a small incision.

Fast, Efficient Procedure

Ready-to-use disposable instruments for a standard core decompression.


Advanced Surgical Technique - Core Decompression of the Femoral Head

Proper surgical procedures and techniques are the responsibility of the medical professional. The following guidelines are furnished for information purposes only. Each surgeon must evaluate the appropriateness of the procedures based on his or her personal medical training and experience. Prior to use of the system, the surgeon should refer to the product package insert for complete warnings, precautions, indications, contraindications and adverse effects. Package inserts are also available by contacting Wright Medical Technology, Inc.

Remove all plastic caps on instruments prior to use.

Step 1 | ACCESS LESION

Use a 2cm stab incision for access. Under fluoroscopic guidance (both AP & LAT views), introduce the 3.2mm fluted guidewire into the lesion. Most lesions are anterior and superior.




Step 2

Introduce the tissue protector over the guidewire and down to the bone prior to drilling.




Step 3 | DECOMPRESS FEMORAL HEAD

Using the 9mm cannulated drill bit, decompress the femoral head by drilling a core approximately 5mm from the endosteal surface of the femoral head. AP and LAT fluoroscopic views should be used to confirm direction.



Step 4 | PLACE WORKING CANNULA

Maintain placement of the tissue protector and remove the drill bit and guidewire. Place the working cannula with obturator through the tissue protector and into the core. Position the working cannula up into the core several centimeters (fit should be snug). Remove the tissue protector and obturator.



For Standard Debridement
with PRO-DENSE® Core Decompression Procedure Kit

Step 5 | DEBRIDE DEAD BONE

Standard debridement can be accomplished using the curette and/or the fluted guidewire.



Note: If not using the X-REAM® tool, go to step 9.

Advanced Core Decompression X-REAM® Percutaneous Expandable Reamer

Complete Steps 1 - 5 with the PRO-DENSE® Core Decompression Procedure Kit Instruments prior to debridement.


Step 6 | ADVANCED DEBRIDEMENT

Advanced debridement can be carried out using the X-REAM® Percutaneous Expandable Reamer.




REAMER PLACEMENT

Introduce the Reamer through the working cannula confirming placement with fluoro (AP & LAT views).



Step 7 | DEBRIDE DEAD BONE



A) Turn the blade control knob $\frac{1}{4}$ turn (or less) clockwise.

Note: It is **extremely important** not to open the blades too far prior to rotating the tool. Otherwise, blade failure will occur.


B) Rotate the entire instrument two full revolutions.

C) Repeat steps A & B until desired expansion is achieved. Use fluoro as needed to monitor the blade expansion.

D) If the blades are opened too far prior to cutting, turn the blade control knob counterclockwise until it stops and withdraw the instrument into the working cannula to collapse the blades. Reinsert the instrument and begin again at Step 7A.

Note: The instrument **must** be withdrawn back into the working cannula in order to collapse the blades.

X-REAM® Tool *in situ* Expansion




WARNING: During expansion, frequently confirm blade position under fluoro (BOTH AP and LAT views). Rotate instrument so blade width can be clearly determined (i.e. blades are perpendicular to view).

CAUTION: Be sure not to violate the subchondral plate during the blade expansion.

Step 8 | REMOVE X-REAM® TOOL

Once debridement is complete, turn the blade control knob counterclockwise until it stops.


Simply withdraw the X-REAM® tool through the working cannula. The blades will self-collapse.



Step 9 | ASPIRATE CORE (RECOMMENDED)

Once debridement is complete, use the suction tip from the PRO-DENSE® Core Decompression Procedure Kit to remove the debrided tissue.

Flushing with a combination of irrigation and suction works best.



Step 10 | GRAFT CORE

Prepare graft per instructions provided in the kit.


Backfill the core using the PRO-DENSE® Injectable Graft (included in kit) to completely fill the surgically-created bone defect.

Begin by placing the needle at the back of the defect and injecting with thumb pressure.

Slowly inject while simultaneously withdrawing needle.

Periodically check graft placement with fluoro.

Slowly remove the working cannula while backfilling the core.



Step 11

Confirm final placement of graft under fluoroscopic guidance and close in standard fashion.



Ordering Information



X-REAM® Percutaneous Expandable Reamer

- 1000-KIT2 X-REAM® Instrument Kit
- 20BL-1200 X-REAM® Blade (Sterile Packed)



PRO-DENSE® Injectible Regenerative Graft

- 87SR-CK15 PRO-DENSE® Core Decompression Procedure Kit 15cc



- 87SR-0410 PRO-DENSE® Injectible Regenerative Graft 10cc
- 87SR-0420 PRO-DENSE® Injectible Regenerative Graft 20cc

REFERENCES

1. SooHoo NF, et al. Cost-effectiveness analysis of core decompression. *J Arthroplasty* 2006; 21(5):670-681.



Wright Medical Technology, Inc.
 1023 Cherry Road
 Memphis, TN 38117
 800 238 7117
 901 867 9971
 www.wmt.com

Wright Medical EMEA
 Atlas Arena, Australia Building
 Hoogoorddreef 7
 1101 BA Amsterdam
 the Netherlands
 011 31 20 565 9060