


ORTHOLOC® 3Di

Ankle Fracture System






SYSTEM OVERVIEW



3Di

 **WRIGHT**™


ORTHOLOC® 3Di: Ankle Fracture System



The Comprehensive Single Tray Solution


Advanced • Versatile • Complete

ORTHOLOC® 3Di: Advanced Design



ORTHOLOC® 3Di Locking Technology

- 30° conical polyaxial locking capability
- Multi-point screw/plate engagement
- Uncompromised polyaxial locking stiffness



Bullet tip plates for submuscular plate insertion

Universal locking holes accept:


- 2.7mm locking
- 3.5mm locking
- 3.5mm non-locking
- 4.0mm non-locking

Compression slots for 3.5 and 4.0mm non-locking screws

Advanced, pre-contoured, low-profile plate geometries

Anatomic screw trajectories for surgical speed and accuracy


ORTHOLOC® 3Di: Comprehensive Implant Solutions




Lateral Fibula



Medial Malleolar




Posterior Fibula



Syndesmosis


Posterior Tibia



1/3 Tubular Plates



Anterior Tibia



Medial Tibia



Anterior Tibia



Anterior Lateral Tibia

O.R. Efficiency

Instrument/Implant Color Coding System



3.5 Cortical Full Thread



4.0 Cancellous Partial Thread




4.0 Cancellous Full Thread



4.0 Cortical Full Thread



4.0 Cannulated




3.5 Locking



2.7 Locking



3.5 Low-Profile





Wright Medical Technology, Inc.

1023 Cherry Road
Memphis, TN 38117
800 238 7117
901 867 9971
www.wmt.com

Wright Medical EMEA

Atlas Arena, Australia Building
Hoogoorddreef 5
1101 BA Amsterdam
The Netherlands
011 31 20 545 0100

™Trademarks and ®Registered marks of Wright Medical Technology, Inc.
DARCO® is a licensed trademark of Wright Medical Technology, Inc.
©2013 Wright Medical Technology, Inc. All Rights Reserved.

008919A 11-Nov-2013